

Stratus EEG™ Quick Reference Guide




LOGIN

Upon login attempt, a two-factor authentication code email will be sent to the email address on file.

PLEASE NOTE: The two-factor authentication code email will come from noreply@stratus EEG.com.

No Email Code? Please check the Spam, Junk, or Clutter folders.

HIPAA regulations require two-factor authentication for all new devices and internet browsers used. Select the **Trust This Device** checkbox to remain authenticated for 30 days per device and browser.




To reset your temporary password:

1. Click the **My Settings** tab
2. Click the **Change password** button
3. Enter your temporary password in the **Old Password** free-text field
4. Enter your unique password in the **New Password** and **Confirm Password** free-text fields
5. Click the **Change** button

To receive two-factor authentication codes by text, enter the preferred mobile number (digits only), including "1" + area code, into the **Mobile Number** free-text field in the Stratus EEG **My Settings** tab. (e.g., 18772183122)

PATIENT VISIT LIST




Search/
Filter

Sort

Patient
List

PATIENT EXAM CARDS





To view the exam properties, hover your mouse over the **Patient Card** and click within the blue highlighted area.

REVIEW AND RETRIEVE

The exam workflow state label is displayed on the bottom left-hand corner of the **Exam Box**.

- Click the **Review** icon to open the exam tracing in a new browser tab.
- Exams are **Archived** after 90 days. To access an archived exam, click the **Retrieve** icon.




(Continued on Back)

PATIENT FOLDER

All patient information and exam details are included within the **Patient Folder tabs** that are categorized by **Patient**, **Study**, and **Interpretation**.

PROPERTIES TAB

The **Properties** tab is automatically displayed when a **Patient Folder** is opened. The Properties tab contains patient demographic information and exam details that can be included in reports. To view, edit, or update the patient information and exam details, click the preferred **Exam Box** within the **Patient Folder**, and then click the preferred **Patient Folder tabs**.



Helpful Hint:
The **Report Events** tab was previously the **Episodes** tab.


PATIENT FOLDER CATEGORIES
Categories are broken down into subcategorized tabs.

EXAM ACTIONS
Preview Reports and manage exams using the **Exam Actions** dropdown menu.

The screenshot shows the StratusEEG interface with the Patient Folder open for exam 72FA55. The 'Properties' tab is active, showing patient demographics, study details, and exam information. A green box highlights the 'Exam Actions' dropdown menu, which is expanded to show options like 'Preview Report', 'Delete', 'Merge', 'Reassign', 'Export', 'Archive', 'Detection', 'Delete Video', and 'Send Exam Link'. Another green box highlights the 'Report Events' tab in the left sidebar. A third green box highlights the 'Patient Folder' categories in the left sidebar.

DOCUMENTS TAB

Documents and reports are displayed in the **Documents** tab. *Reports* are added to the **Documents** tab when they are previewed and signed using the **Preview Report** option in the **Exam Actions** menu. *Documents* can be uploaded from a computer or electronic device using the **Add Document** button located at the top left-hand corner of the **Documents** tab.



The screenshot shows the StratusEEG interface with the Patient Folder open for exam 72FA55. The 'Documents' tab is active, showing a list of documents. A blue box highlights the 'Add Document' button. Below the document list, three blue arrows point to the 'Download', 'Sign', and 'View' buttons.