



Trackit Plus Acquisition Software

Trackit Plus software, designed for use with all our hardware, is a simple to use, cost-effective, EEG and video acquisition package that supports CPT Review requirements.



Above: The Trackit Plus lock screen on the Trackit Solo ambulatory base station and Trackit T4A

Don't let its retro looks deceive you

Trackit Plus software was originally designed to support ambulatory EEG recorders such as the Trackit Mk2 and Mk3. However its reliability and simplicity of use have made it the EEG acquisition software of choice for many of our clients. With its robust recovery features and tamper-resistant lock screen, its retro looks and familiar MS Office ribbon, its no surprise that it continues to be a firm favorite among EEG Techs and Service Providers alike.

Trackit Plus Software is compatible with:

Trackit Solo	Trackit T4A	Trackit Mk3	SleepWalker
Trackit T4	R40	Photic	TrackitCam
All PCU options (Trackit and T4)		Incereb Neon Kits	

Highlights



Bluetooth
Wireless



USB
Interface



Ambulatory
Recording



In-Home
Recording



Lock Screen



Recovery
Features



Summary of Features and Benefits

Key Benefits

- Support for all Lifelines hardware
- Familiar Office Ribbon UI
- USB and wireless Bluetooth interface
- Ambulatory recording, PC recording or both
- Built-in User Manuals
- Lockscreen

Data Capture

- 16 bit edf support
- 24 bit bdf support, including DC
- Up to 80 channels
- Up to 2kHz sampling
- >72 hour recording support

Video Capabilities

- Video synchronised recordings
- Mp4 and avi video formats
- Configurable Video Resolution
- Playback/review of edf/bdf files with synchronised video
- Preferred software for the TrackitCam

Settings

- Calibration
- High filter, low filter and notch filter
- Variable display sensitivity
- Variable chart speed
- Variable playback speed
- Individual channel control
- Signal and Recording Setup Editors
- Electrode Hookup display
- Measurement graticule
- SpO2 support, including HR and plethysmograph

Additional Features

- Bluetooth out of range alert
- Auto resume record after power off
- Montage editor, including 10/20 and 10/10 support
- Photic setup editor - manual, sweep and complex patterns
- Hyperventilate
- Electrode impedance measurement, incl. Ref electrode
- Continuous electrode impedance
- File manager utility for assembling files after study completion
- Live Event Placement
- Configurable events, including free-text event