


ACCESSING & LOGGING INTO STRATUS EEG:


Login Webpage Information



New Stratus EEG Users: Your Stratus EEG username and temporary password will be sent to your preferred email address.

Centrum Desktop Application Users converting to Stratus EEG: Your Stratus EEG username and password are the same as your Stratus Centrum & Review Desktop Application username and password and will be emailed to your preferred email address.

PLEASE NOTE: There must be an updated email address associated with your Stratus EEG user account to convert from the Centrum Desktop Application to Stratus EEG Web seamlessly.



Login

Username

Password

Remember My Login

Login

[Forgot Password](#)

Forgot password?

Click the [Forgot Password](#) option to reset your password.

PLEASE NOTE: After 5 unsuccessful login attempts your Stratus EEG user account will be locked. If you need your account unlocked, email support@lifelinesneuro.com.

OPTIONAL - Remember My Login:

Click the **Remember My Login** checkbox to initiate the auto-fill, save password option in Microsoft Edge.