



# MY SETTINGS TAB:

## US Physician Toggle – Physician Information for Reports



The screenshot shows the 'My Settings' page in the StratusEEG application. The left sidebar contains navigation options: Patient Visits, Tasks, My Tags, EEG Settings, and My Settings (selected). The main content area is titled 'My Settings' and includes the following fields and controls:

- User Name: stratus (with a 'Change password' button)
- Full Name\*: Stratus Demo User
- Email Address\*: productmanagement@stratusneuro.com
- Mobile Number: (181)798-0545
- Receive email notifications
- Upload Signature button
- US Physician** toggle (highlighted with a green box)
- NPI Number\* (text input)
- Street Address\* (text input)
- City\* (text input)
- State\* (dropdown menu) and Zip Code\* (text input)
- Credentials (dropdown menu)
- Specialty (dropdown menu)
- Cancel and Save buttons at the bottom.



Click the **US Physician** Toggle in the Stratus EEG, **My Settings** tab to enter Physician details.

**PLEASE NOTE:** This information will be included when using the **Preview Report** features from the **Patient Folder**, **Exam Actions** dropdown menu.