

# Trackit Plus updates



Trackit Plus 2.8.4.7 brings updates that are optimized for our users conducting at-home long term video EEG's.



A new lock screen feature has been implemented, allowing you to lock the software with a password that prevents the software from being exited or closed by unauthorized users. This lock screen is a great tool that still allows quick access to the most commonly used features during a recording, like event placement.



Upon stopping a study, our latest update now includes a convenient built-in file manager that simplifies the process of retrieving amplifier data and merging it with any relevant video files or recorded events saved on your PC. This feature is specifically designed to enhance the efficiency of data management and streamline your workflow.

## Other features include:

- Support for dual cameras
- Amplifier battery percentage logged as event at the beginning of a study
- Compatibility to import BDF+ files for review

Discover the full potential of these new features with our Technical Support Hub, featuring comprehensive how-to-guides and informative videos.

*For more information on updating your Trackit Plus software, please contact our support team via phone or email.*



Imagine EEG Anywhere

support@lifelinesneuro.com  
866-889-6505

lifelinesneuro.com

A **SIDIS** Company